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ENGLISH

CONTENTS

I - INTRODUCTION	6
II - CAUTIONS	6
III - DESCRIPTION	6
3.1 ITEMS DELIVERED WITH THE PERIO KIT:	6
3.2 OVERVIEW	6
3.3 TECHNICAL CHARACTERISTICS:	7
IV - INSTALLATION / SWITCHING ON	7
4.1 UNPACKING the device	7
4.2 INSTALLATION	7
4.3 FIRST USE	8
V - AIR-N-GO PERIO USAGE PROTOCOL	8
5.1 USING THE AIR-N-GO HANDPIECE	8
5.2 USING THE PERIO NOZZLE	10
VI - SERVICING AND STERILIZATION	12
6.1 RECOMMENDATIONS:	12
6.2 SERVICING AND CLEANING OF THE PERIO NOZZLE	13
6.3 SERVICING AND CLEANING THE TANK	15
6.4 SERVICING AND CLEANING THE HANDPIECE	18
VII - CHECKING / MAINTENANCE	18
7.1 CHECKING	18

7.2 MAINTAINING THE “STOP POWDER” PISTON (FIG.5-8)	19
7.3 MAINTAINING THE METAL SHAFT	19
7.4 OPERATING FAULTS	19
VIII - MANUFACTURER RESPONSIBILITY	19
IX - REGULATIONS	20
X - MANUFACTURER	20
XI - DATE OF REVISIONS OF THE USER MANUAL:	20
XII - FIGURE	21

I - INTRODUCTION

Created by SATELEC®, the AIR-N-GO™ PERIO kit is compatible exclusively with the AIR-N-GO device, transforming it into a dental air-polishing device used to carry out the following types of treatment:

- subgingival treatments of periodontitis or peri-implantitis, exclusively using the glycine-based 3M ESPE® Clinpro™ Prophy Powder¹.

To maximize use of this cutting-edge product for as long as possible, please read this user manual carefully before switching on, using and servicing the device.

II - CAUTIONS

CAUTION: The United States Federal Law restricts the use of this device to qualified dental health professionals only or professionals capable and qualified or under their control. To reduce the risk of accidents, please read and adhere to the following precautions:

Device operators:

- Use of the AIR-N-GO device is limited to qualified dental health professionals only, capable and qualified within the usual scope of their activities.
- If you have received this device by mistake, contact its supplier to arrange for its return.

Precautions of use:

- The AIR-N-GO PERIO kit is an option exclusively associated with the AIR-N-GO device. It is designed to be exclusively used with the 3M ESPE glycine-based Clinpro™ Prophy Powder.

III - DESCRIPTION

3.1 ITEMS DELIVERED WITH THE PERIO KIT:

- A graduated PERIO nozzle that can be adjusted and sterilized (Figure 1)
The PERIO nozzle comes with 3-6-9 mm depth markings that allows rapid assessment of the health of a patient during periodontitis or peri-implantitis subgingival treatments.
- A green PERIO “CLIP-N-GO” powder tank (Figure 2)
- A PERIO powder tank lid (Figure 2)
- A 100g bottle of 3M ESPE Clinpro™ Prophy Powder (Figure 3)
- A PERIO nozzle cleaning probe (Figure 4)

COMMENT:

To benefit from the AIR-N-GO polishing device PERIO option, assemble the PERIO tank and nozzle delivered with your AIR-N-GO PERIO kit.

The AIR-N-GO handpiece is required to correctly operate these parts.

3.2 OVERVIEW

Using PERIO, the AIR-N-GO handpiece comprises:

- A graduated PERIO nozzle that can be adjusted and sterilized (Fig.5-1)

¹ Clinpro™ Prophy Powder is a patented trademark of 3M ESPE.

- A metal head (Fig.5-2)
- A handpiece body comprising:
 - A non-slip silicone ring (Fig.5-3)
 - A plastic front body (Fig.5-4)
 - A rear body (Fig.5-5)
 - A turbine adapter depending on the dental chair model (Fig.5-9)
- A green PERIO “CLIP-N-GO” powder tank (Fig.5-6)
- A PERIO powder tank lid (Fig.5-7)

The maximum powder level is marked by a “MAX” line (FIG. 6).
The “STOP POWDER” piston used to empty the system is located on the tank lid (Fig.5-8).

3.3 TECHNICAL CHARACTERISTICS:

1. Device

 Refer to the AIR-N-GO device user manual.

2. Powders

 Refer to the AIR-N-GO powder instructions provided in the AIR-N-GO box.

IV - INSTALLATION / SWITCHING ON

4. 1 UNPACKING THE DEVICE

Upon receipt of the device, check for any damage that may have occurred during transportation and contact your supplier if necessary.


4.2 INSTALLATION

 **Important:**

- For the subgingival treatment of periodontitis and peri-implantitis: only use the PERIO nozzle with 3M ESPE glycine-based Clinpro™ Prophy Powder and the green tank.
- Make sure there are no signs of moisture around the nozzle or in the powder tank. Eliminate any moisture found by wiping and drying with the multipurpose syringe air spray.

To transform the AIR-N-GO device, based on PERIO, you need to:


- Assemble the PERIO nozzle on the AIR-N-GO handpiece body, manually screwing on the nozzle clockwise.

- Make sure that the nozzle is securely in place.


- Install the green tank on the handpiece body (Cf. 6.3 b) in accordance with the tank assembly procedure chap.VI.

-Connect the device by simply plugging it into the turbine cord connector.

4.3 FIRST USE

 **The tank and its lid cannot be sterilized.**

Recommendations:

The use of the AIR-N-GO device is limited to qualified dental health care professionals, capable and qualified within the usual scope of their activities (dental surgeons, hygienists, etc.).

Prior to first use, sterilize the nozzle in an autoclave in accordance with these parameters:

- Autoclave: Type B in compliance with standard EN 13060.
- Sterilization temperature: 134° C.
- Sterilization duration: 18 minutes.
- Pressure: at least 2 bars.

The user instructions for the device must be complied with (Refer to Chap V).

V - AIR-N-GO PERIO USAGE PROTOCOL

5.1 USING THE AIR-N-GO HANDPIECE

AIR-N-GO can be implemented as follows:


- Remove the bayonet system tank from the handpiece by rotating it a quarter of a turn anti-clockwise.

- Unscrew and remove the lid from the powder tank.

- Fill the tank with the powder recommended by SATELEC making sure that you do not fill it above the "MAX" line.


- Put the lid back on the tank.

- Put the tank back on the handpiece by rotating it quarter of a turn clockwise.

- Connect the handpiece to the turbine cord connector.


- Adjust the water flow (15 ml / min \pm 5 ml / min).
- Press the control pedal on the dental chair.
- Treat your patient as required.


⚠ At the end of treatment, or prior to sterilization, simultaneously press the "Stop Powder" piston on the tank lid (Fig.1-8) and the control pedal on the dental chair to bleed the air system.

- Remove the tank from the handpiece by rotating it quarter of a turn counter-clockwise.

- Unscrew and remove the lid from the powder tank.

- Empty all the powder from the tank into the chair spittoon or the sterilization room sink.


- After each use, it is important to service and/or sterilize all equipment in accordance with the procedures defined in chapter VI.

5.2 USING THE PERIO NOZZLE

 The patient and dentist must wear suitable protection accessories (glasses, masks).
- Apply Vaseline® to the patient's lips before polishing.

- Screw the PERIO nozzle onto the AIR-N-GO device.


- Fill the PERIO tank with 3M ESPE Clinpro™ ProphyPowder only.

- Next, assemble the PERIO tank by rotating it quarter of a turn clockwise.

- Connect the AIR-N-GO PERIO device by simply plugging it into the turbine cord connector.


⚠ The AIR-N-GO PERIO kit is designed exclusively for the subgingival treatment of periodontitis and peri-implantitis.

- Use a surgical suction cannula (as powerful as possible) and hold it near the area to be treated (reduces bacterial air contamination).

- Treat your patient as required. Each tooth or implant must be treated individually.

- Insert the nozzle to the bottom of the pocket exerting the same amount of pressure as with a periodontal probe.

- Activate AIR-N-GO by pressing the control pedal, only when the PERIO nozzle is in the periodontal pocket.

- Gently move the nozzle upwards moving it slowly from side to side without pressure for 5 seconds.

- Be careful not to damage the epithelial attachment.

- Leave AIR-N-GO™ on even when you have gone past the top of the pocket. The gingival festoon must be gently treated in order to obtain a substantial reduction of the pocket.

- The air/powder spray continues for a few seconds after releasing the control pedal. Take this into account and wait before removing the nozzle from the patient's mouth to prevent irritation of the mucosa.

- Repeat the procedure on the other surfaces of the tooth or implant, if necessary, or for all pockets and surfaces to be treated.

VI - SERVICING AND STERILIZATION


6.1 RECOMMENDATIONS:

a) At the end of the working day, you must:

- Simultaneously press the “STOP POWDER” piston on the tank lid (Fig.1-8) and the control pedal on the dental chair to bleed the air system.

- Disconnect the handpiece from the turbine connector.

- Completely empty the powder tank and clean it with compressed air.

- Unscrew and clean the nozzle.

b) After each use of the device or its accessories:

- The device must be disconnected during cleaning/disinfection procedures.

⚠ The tank and its lid cannot be sterilized.

- Before starting the cleaning procedure, it is essential to disassemble the tank and its lid. See 6.3.a

- If the device has not been used or has been misplaced for a period of time, you will need to empty all powder from the tank into the chair spittoon or the sterilization room sink and clean the tank.


- Do not use the device if it appears to be faulty.
- Sterilized parts must be left to cool and dry until they reach ambient temperatures prior to being used again.
- Before sterilization, check the cleanliness of your autoclave and the quality of the water used.
- After each sterilization cycle, immediately remove the parts from the autoclave to reduce the risk of corrosion of metal parts.
- Do not use an abrasive product to clean the device.

6.2 SERVICING AND CLEANING OF THE PERIO NOZZLE

! Important:

During sterilization, metal parts of a different nature must not be in contact with each other. Any contact will lead to the creation of galvanic couples resulting in local deterioration. To prevent this phenomenon, place the objects individually in a sterilisable microporous bag or cloth.


a) Preliminary disinfection / Chemical cleaning

- After each treatment, unscrew the PERIO nozzle from the body of the AIR-N-GO handpiece.

 Remove any residual powder from the PERIO nozzle using the PERIO nozzle cleaning probe provided with your AIR-N-GO PERIO Kit.

- Only use the cleaning probe supplied in your AIR-N-GO PERIO Kit.


- Remove any residue from inside the PERIO nozzle through the front and the rear using the air spray of the chair multipurpose syringe.

- Unscrew the head from the AIR-N-GO handpiece body.

- Immerse the PERIO nozzle in an ultrasonic tank containing Listerine®² for at least 10 minutes.

- Remove the nozzle and tap it on a hard surface to remove any remaining particles.
- Rinse the PERIO nozzle under running water for at least 30 seconds.

² Listerine® is a patented trademark of Pfizer

b) Drying

- Using the multipurpose syringe air spray and/or a clean disposable non-woven cloth, dry until there is no liquid left.

c) Sterilization

CAUTION:

- Prior to sterilization, remove any liquid from the accessories.
- To ensure sterile and aseptic conditions, keep all parts in hermitically sealed containers or bags designed for medical use.
- After each treatment, store the disassembled nozzle separately in disposable sterilization sleeves or bags, in compliance with the specifications defined in standard EN ISO 11607-1 or in any equivalent standard that may be required by national regulations.

- Sterilise the nozzle in compliance with the following parameters:
 - Autoclave: Type B in compliance with standard EN 13060.
 - Sterilization temperature: 134°C.
 - Sterilization duration: 18 minutes.
 - Pressure: 2 bars minimum

The head must be removed from the handpiece prior to any sterilization procedure.

d) Storage

- Next, store the sterilized products in a dry place, free from dust and at ambient temperature.
- Before the next use, in the event of a packaging non-conformity, repackage and resterilise in accordance with the defined protocol.
- In the event of visible contaminations in the bag, destroy the product.

e) Product disposal


Dispose of the nozzle in bins reserved for waste following treatment with risks of infection.

6.3 SERVICING AND CLEANING THE TANK


a) Tank disassembly procedure

- Firmly hold the AIR-N-GO handpiece in one hand.
- Hold the tank (Fig.1-6) with the other hand.
- Rotate the tank one quarter of a turn anti-clockwise to release it.

- Pull it to release it from its recess.

b) Tank assembly procedure

! Prior to each tank installation, we strongly advise you to remove all traces of powder using the multipurpose syringe air spray or disinfectant/cleaning wipes containing alcohol, amphoteric and biguanide.

- Fill the PERIO tank with 3M ESPE Clinpro™ Prophy Powder only.
- Reposition the tank in its recess by pushing it.

- Rotate the tank a quarter of a turn clockwise to secure it by locking the bayonet system.

- Check that it is securely in place.

c) Cleaning and maintenance of the bayonet tank:

We recommend that you clean the tank after each use.

- Fully and systematically empty all powder from the tank if the device is not going to be used for several hours (e.g. overnight). Residual air humidity alters the properties of the powder.

- Inside of the tank and the lid:

- Remove all traces of powder residue from inside the tank and the lid using the multipurpose syringe air spray.

- Clean the tank openings using the PERIO probe (Fig.4).

- Clean the tank's central outlet openings using the probe (Fig.7)


- Fully eliminate all traces of powder residue by using the multipurpose syringe air spray again and placing your index finger over the tank outlet.

- Use alcohol-based cleaning/disinfectant wipes to wipe the inside of the tank and its lid.

- Completely dry the powder tank and its lid using the multipurpose syringe air spray.

! The presence of moisture may cause blockage of the device powder or air circulation pipe.


- Use alcohol-based disinfectant/cleaning wipes to clean the inside of the tank and its lid.
- The inside of the tank must be dried using the multipurpose syringe air spray.
- Make sure that the powder tank is completely dry.

- Outside of the tank:

- The outside of the tank and its lid must be cleaned and disinfected daily with disinfection products or disinfectant/cleaning wipes containing alcohol, amphoteric and biguanide.

6.4 SERVICING AND CLEANING THE HANDPIECE

Refer to the AIR-N-GO device user manual for advice on the servicing and cleaning of the AIR-N-GO handpiece.

VII - CHECKING / MAINTENANCE

! Important:

In the event of a fault, we recommend that you contact the supplier of your device rather than a general repairer who could render your device dangerous for you and your patients.

7.1 CHECKING

You will need to check the device and its accessories regularly to detect any fault or deterioration. Where necessary, replace faulty parts.

 **CAUTION:** you must comply with this advice, notably to prevent any risk of clogging (by the powder) inside the AIR-N-GO handpiece.

7.2 MAINTAINING THE “STOP POWDER” PISTON (FIG.5-8)

- Make sure that the “Stop Powder” piston on the tank lid is working correctly (Fig.1-7) by pressing and releasing it a few times.

- Check that there is an o-ring in the lid.

- Check the status of wear of the o-ring and replace it if necessary.

- Spare o-ring kit to order: F10205

7.3 MAINTAINING THE METAL SHAFT

For details, refer to the AIR-N-GO device user manual.

7.4 OPERATING FAULTS

For details, refer to the AIR-N-GO device user manual.

VIII - MANUFACTURER RESPONSIBILITY

The manufacturer shall not be held responsible in the event:

- of non-compliance with manufacturer recommendations during installation.
- of intervention, modifications, or repairs carried out by people who have not been authorised by the manufacturer.
- of usage other than that specified in this manual.
- of the use of accessories and/or powders that have not been recommended by SATELEC.
- of the non-compliance of usage conditions.

Note: the manufacturer reserves the right to modify the device and/or the user manual without notice.

IX - REGULATIONS

This medical device is classified IIa in accordance with the applicable European directive pertaining to Medical Devices.

This equipment was designed and manufactured in accordance with an EN ISO 13485-certified quality assurance system.

X - MANUFACTURER

SATELEC

17 AVENUE GUSTAVE EIFFEL - BP30216

33708 MERIGNAC CEDEX - FRANCE

Tel: +33 (0) 556.34.06.07

Fax: +33 (0) 556.34.92.92

E.mail: Satelec@actongroup.com

www.acteongroup.com

	Refer to the accompanying documentation
	Sterilization at 134°C in an autoclave
	EC marking
AAAA	Year of manufacture

XI - DATE OF REVISIONS OF THE USER MANUAL:

This user manual was last revised in 10/2010.

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