




Technical Datasheet

Désignation	Lubricant Spray R&S				
Référence	<table><tr><td>Lubricant Spray</td></tr><tr><td>3149969</td></tr></table>			Lubricant Spray	3149969
Lubricant Spray					
3149969					
Product description	Lubricant for rotative instruments				
Product picture					
Composition	HYDROCARBON PROPELLANT REFINED MINERAL OIL PENTANE				
Packaging	<ul style="list-style-type: none">- 500 ml Spray- Container under pressure.				
Storage	<ul style="list-style-type: none">- Harmful to aquatic life with long lasting effects.- Protect from heat above 50°C- Do not puncture or burn can- Keep away from sources of ignition.- No smoking. Do not breathe vapour/spray. Use only in well ventilated areas.- This material and its container must be disposed of in a safe way.- Keep out of the reach of children.				

Version 3.0

Mis à jour 2015/12/02 par Laetitia Decaudin

CFPM - 2 bis chemin du Loup 93290 Tremblay-en-France - FRANCE - Tel. +33 (0)1 45 91 31 82