



A PARTIR DE:

**5€
,10**

Bombillas para aparatos de iluminación quirúrgica DEGRÉ K

DESIGNACIÓN	REF.	EUROS
°K201	15-848	41,00
°K203	15-850	11,10
°K204	15-851	21,70
°K205	15-852	40,70
°K206	15-853	7,60
°K208	15-855	6,15
°K209	15-856	5,10
°K211	15-858	12,50
°K212	15-859	5,50
°K214	15-837	44,00
°K215	15-834	42,70
°K217	15-835	9,80
°K220	15-619	20,90
°K222	15-617	67,00
°K230	15-860	5,40
°K233	15-441	19,20
°K234	8-044	20,60

MARCA	MODELO	VOLTIOS	VATIOS	N° °K
ADEC	6300	24	150	°K209
	Performer	12	55	°K208
ANTHOS	Venus	17	95	°K233
BELMONT	24 V	12	55	°K208
		17	70	°K221
CASTELLINI	Vega, Monta, Solar	12	50	°K212
	LUNA	12	75	°K217
EURODENT	Isolite (300 horas)	17	150	°K203
	Isosun	17	100	°K211
	Riunito	17	150	°K210
	Isolite (400 horas)	17	150	°K209
	70, B4843, Sunlight	17	150	°K203
	S2000, S70, S90	17	95	°K204
FARO	70, B4843, Sunlight	17	150	°K209
	EDI	17	95	°K233
	Forest, Koolbeam, S2000, S70, S90	17	150	°K234
FEDESA	Lineo, Midway, Scorpio	17	95	°K204
FLEX		12	55	°K208
	1400	17	150	°K203
KAVO	1400	17	150	°K209
	1400, 1410A, Esttica 1040	17	150	°K214
	Vision	17	50	°K230
PLANMECA	Delight	24	150	°K209
	Delight	24	150	°K203
	PM 2002	24	70	°K221
RITTER	D44 (casquillo, 2 pasadores) - D 428	24	150	°K210
	Ritter 130, Lumastar new	24	150	°K209
	Star track, "K" Super starlite	24	150	°K234
	Sirolux S 24 V, Fantastic (C1,C2,C3,C4)	24	150	°K220
SIRONA	Sirolux S 12 V	12	100	°K206
	Sirolux	12	55	°K208
	Sirolux E (C6,C8)	12	50	°K230
VITALI	Vitali Cammbologna	17	95	°K204
XO	XO4	17	50	°K230

El No 1 de los MEJORES PRECIOS