

<b>Designation</b>	<b>Saliva ejectors – removable tips</b>					
<b>Reference(s)</b>	<b>Colours</b>	Translucent	Blue	Green	Orange	Pink
	<b>Reference</b>	3127251	3127252	3127253	3127254	3127255
<b>Product picture</b>						
<b>Product description</b>	<b>Formatting easy</b>	because armed with a copper metallic wire				
	<b>Tips</b>	removable, little, smooth, rounded				
	<b>Use</b>	ideal for all uses, even on irritated mucous membranes				
	<b>Diametre</b>	6.5 mm				
	<b>Length</b>	15 cm				
<b>Characteristics and compositon</b>	Made from PVC					
<b>Packaging</b>	<ul style="list-style-type: none"> <li>- Store at room temperature in a dry and ventilated area</li> <li>- Plastic bag of 100 units</li> </ul>					
<b>Shelf life</b>	n.a					
<b>Indications of use</b>						
<b>Legal information</b>	Medical Device of class IIA under European Directive 93/42/EEC including Directive 2007/47/CE CE 0051 Marking- Organismo de certificación: IMQ S.p.A Manufacturer: EURONDA					

<b>Designation</b>	<b>Saliva ejectors – fixed tips</b>		
<b>References</b>	<b>Colours</b>	transparent	white
	<b>Reference</b>	3135200	3135201
<b>Product Picture</b>			
<b>Product description</b>	<p>Easy shaping because of copper wire reinforcement</p> <ul style="list-style-type: none"> <li>- Tips : small, smooth and rounded</li> <li>- Ideal for all uses, even on irritated walls</li> <li>- Diameter : 6,5 mm</li> <li>- Length : 15 cm</li> </ul>		
<b>Characteristics and composition</b>	Made from PVC		
<b>Packaging</b>	<ul style="list-style-type: none"> <li>- Store at room temperature in a dry and ventilated area</li> <li>- Plastic bag of 100 units</li> </ul>		
<b>Shelf life</b>	5 years from the date of manufacture		
<b>Indications of use</b>			
<b>Legal information</b>	<p>Medical Device of class IIA under European Directive 93/42/EEC including Directive 2007/47/CE            CE marking- 0051-IMQ S.p.A            Manufacturer: ASA DENTAL</p>		