

1187 SAFETOUCH® ADVANCED BLACK NITRILE POWDER-FREE Class I, Medical Device and Category II, Personal Protective Equipment

FEATURES & BENEFITS

- Excellent grip objects thanks to its elasticity and dexterity
- Its textured fingers improve grip in dry and wet conditions.
- Watertight and airtight

Recommendation:

- Glove for professional use
- In the healthcare and industrial environment, this glove is used to make barrier against microorganisms between the wearer of the glove and the examined people or objects / materials handled.
- Powder-free, it reduces cross contamination and limits the risks of irritation.
- For food handling / preparations (HACCP)



PRODUCT DESCRIPTION

Product name :	SafeTouch® Advanced Black
Product type:	Single use, non-sterile
Material :	Nitrile
Colour :	Black
Shape:	Ambidextrous
Cuff:	Rolled
Surface:	Textured fingertips
Powdered :	No
Length:	240 mm
Palm thickness :	0,1 ± 0,03 mm
Finger thickness :	0,12 ± 0,03 mm
Weight:	5,2 g
Country of origin:	Malaysia
AQL :	< 1,5
Tensile Strength :	
✓ Before aging	≥14MPa
✓ After aging	≥14MPa
Elongation:	
✓ Before aging	≥ 500%
✓ After aging	≥ 400%

Item	Size
1187-A	T 6 XS
1187-B	T 7 S
1187-C	T 8 M
1187-D	T 9 L
1187-E	T 10 XL



INDUSTRIES

Food / Non residential catering / Wellness / Security / Community / Administration / Health / Medical / Sanitary / Maintenance / Cleaning / Cleaning room / Lab / Pharmaceutical / Research

1187 SAFETOUCH® ADVANCED BLACK NITRILE POWDER-FREE

Class I, Medical Device and Category II, Personal Protective Equipment

CERTIFICATION

In accordance with requirements of medical devices Dir.93/42/EEC modified by Dir.2007/47/EC
Comply with the EN 455-1/3/4 and ASTM D6319 for single use medical gloves

Manufactured under ISO 13485 and ISO 9001 certified quality system

In accordance with requirements of personal protective equipment Dir.89/686/EEC

Comply with the EN420 & EN 374-2: 2003 for protection against microorganisms

In accordance with requirements of Reg.(EC) 1935/2004 and French decree of November 9th, 1994.

This glove is suitable for food contact.

*Restriction : Contact with all types of food excluding acidic food stuffs (pH<4,5) and fatty food stuffs (see Annex 1).

STORAGE CONDITIONS

Normal conditions of conservation & storage: not to be exposed to moisture and sunlight, must be stored at a temperature between 5°C and 40°C

Shelf life of the product: 5 years

LOGISTIC INFORMATION



Specifications outer case



Specifications inner case

Colour code dispenser and carton

Item	Size mm	Gross weight kg	QTY/ Pallet
1187-A	340 x 250 x 250	5	63
1187-B	340 x 250 x 250	5,8	63
1187-C	340 x 250 x 250	6,1	63
1187-D	340 x 250 x 250	6,9	63
1187-E	340 x 250 x 250	7,2	63

Size mm	QTY
240 x 120 x 65	10 x 100
240 x 120 x 65	10 x 100
240 x 120 x 65	10 x 100
240 x 120 x 65	10 x 100
240 x 120 x 65	10 x 100

Colour

EAN Codes

Article	Boîte distributrice	Carton
1187-A	3 700 381 514 758	3 700 381 514 369
1187-B	3 700 381 514 765	3 700 381 514 376
1187-C	3 700 381 514 772	3 700 381 514 383
1187-D	3 700 381 514 789	3 700 381 514 390
1187-E	3 700 381 514 796	3 700 381 514 406