



TECHNICAL SHEET

Product Name	Nerve broaches R&S				
Reference	Nerve broaches N°	Colors R&S	Colors ISO	Available length	Reference
	15	●		21 mm	0281877
	20	●	●		0281878
	25	●	○		0281879
	30	●	●		0281880
	35	○	●		0281881
	40	●	●		0281882
	45	●			0281883
	50	●	○		0281884
60	●	●	0281885		
Visual					
Description	The nerve broaches allow you to remove the pulp.				
Indications	Single use for all numbers				
Physical properties	<p><u>Geometrical specifications:</u> Nerve broaches are instruments, with a conicity of 1% and barbed around the circumference. The instruments comply with the applicable paragraphs of ISO 3630-4.</p> <p><u>Materials used:</u> -Handle: Metal -Blade: Stainless Steel</p> <p><u>Pre-disinfection, cleaning and sterilization:</u> Instruments sold in non-sterile packaging must be pre-disinfected, cleaned, and sterilized at 134°C for 18 minutes before use.</p> <p>Nerve broaches should not be immersed in soda products.</p>				

Packaging / stockage	- Sold in packs of 12 non-sterile instruments for each number.
EC Class	Class IIa medical device. CE 1639 marking. Certifying body SGS BELGIUM. Manufacturer : CFPM. Code GMDN 40529

Version 3.0
Updated 22/10/2020 by Erine TCHAKOUNTE