


# « Ultralight » Endo Holders

**New !**


**SPECIALLY DESIGNED FOR WASHER DISINFECTOR USE**



***New Concept, New Design !***

*Very light, those Ultralight Endo Holders have been created specifically for use in washer disinfectors, or in ultrasonic bath.*

- When it is closed, the Ultralight endo holder is a cassette, ensuring a high protection of your endodontic instruments.
- Shiny finish (electro-polishing), avoiding brown marks
- 3 models are available :
  - **Ultralight Endo 4** : for 4 endo files
  - **Ultralight Endo 8** : for 8 endo files
  - **Ultralight Endo 12** : for 12 endo files



- The internal holder can be rotated : when opening the cassette, you turn the internal holder vertically to present the products. This allows a better handling of the instruments.
- The small silicone tubes allow a secure fitting of your instruments.
- Le sterimeter allows you to count the number of uses of the sequence.

Made of Stainless Steel and silicone, those holders have been designed for use with all disinfection and cleaning solutions, as well as all sterilization processes.

REF	DESIGNATION		SIZE
182320	Ultralight Endo 4	4 instruments	50 x 59 x 19 mm
182380	Ultralight Endo 8	8 instruments	84 x 59 x 19 mm
182390	Ultralight Endo 12	12 instruments	119 x 59 x 19 mm

