



Designation	Sterilisation pouches Medibase									
Reference(s)	SELF SEALING POUCHES									
	Size	8,5 x 16,5cm	7 x 25,5cm	9 x 25,5cm	14 x 28cm	25,5 x 38cm	5,7 x 13cm	18 x 33cm	9 x 13cm	30 x 45 cm
	Ref	3135469	3135470	3135471	3135472	3135473	3135474	3135475	3135476	3135477
	HEAT SEALING POUCHES									
	Size	7X25,5CM			9X25,5CM			20X33 AVEC SOUFFLET DE 5,5CM		
Ref	3135479			3135480			3135478			
Product picture										
Product description	Sterilization bags recommended for autoclaving or chemiclave.									
Characteristics and composition	Paper	The basis paper weight is $\pm 60\text{g/m}^2$, fitting the requirements of sterile packaging manufactured paper which list on ISO11607-1, 2 and EN868-5; According to the pouch designs, the printed inscription on the rolls is navy blue (Pantone298U), printing steam and EO indicator.								
	Film	Lucid wathet blue hyperthermal PET/CPP complex film, thickness is $\pm 0,052\text{mm}$, weight is 56 g/m^2 , conforming to the related requirements of EN868-5.								
	Ink	Using ink; the min printing area of indicator should be 100mm^2 , may not fade or become stained.								
Packaging	Box of 200 pouches, available in 10 sizes									
Shelf life	5 years from manufacturing date									
Indications of use										
Legal information	Medical Device of class I under European Directive 93/42/EEC including Directive 2007/47/CE CE Marking Manufacturer : FOMED Industries Inc									