



Filtek™

Dental Restoratives

Experience the difference



That was then.

Is it time to update from the composite you're using now?

3M™ Z100 MP Restorative



- Microhybrid to enable universal use
- Developed to overcome the esthetic challenges of traditional macrofills
- Overcomes the wear and strength challenges of microfills


3M™ Filtek™ Z250 Universal Restorative



- Microhybrid to enable universal use
- Improved resin system resulting in better handling
- Available in 12 shades

In their day, microhybrids offered a significant improvement over older macrofills and microfills. Since then, 3M has raised the bar with our TRUE nanotechnology.

Microhybrid technology



As the resin wears, large and unevenly sized filler particles are exposed. This creates a rough and uneven surface as they break away from the resin matrix.

What does this mean for your restoration?

- Compromised esthetics
- Decreased polish retention

This is now.

3M™ Filtek™ Supreme XTE Universal Restorative



Exceptional esthetics and outstanding strength

- TRUE nanofiller technology enables unsurpassed esthetics and polish retention
- Over 600 million restorations globally*
- Outstanding strength and esthetics for both anterior and posterior use
- Exceptional handling
- Excellent wear resistance
- Available in 36 shades, 4 opacities for exceptional esthetics

3M™ Filtek™ One Bulk Fill Restorative




Simple, one-step bulk fill placement up to 5 mm

- Two innovative monomers enable bulk filling by providing stress relief
- TRUE nanofiller technology for excellent strength and wear resistance
- Fast and easy posterior restorations
- One-step placement, no additional capping layer
- Increased opacity for improved esthetics
- Excellent handling

The advancement of our TRUE nanofiller technology has resulted in a full line of restorative solutions that offer unsurpassed esthetics and simplicity.

TRUE nanofiller technology



Nanoparticles are uniquely formed into nanoclusters. As these nanoclusters wear, they maintain a uniform and smooth surface.

What does this mean for your restoration?

- Excellent esthetics
- Higher polish
- Improved polish retention
- Excellent wear resistance

*3M™ Filtek™ Supreme Universal Restorative line.

Ordering Information

3M™ Filtek™ Supreme XTE Universal Restorative



Item #	Product Information
4911S	3M™ Filtek™ Supreme XTE Universal Restorative Single Shade Intro Kit—Capsule Kit includes: 40 – 0.2g capsules; 20 each of Body shades A2, A3; Technique Guide; Instructions for Use
4910S	3M™ Filtek™ Supreme XTE Universal Restorative Single Shade Intro Kit—Syringe Kit includes: 4 – 3g syringes; 2 each of Body shades A2, A3; Technique Guide; Instructions for Use

3M™ Filtek™ One Bulk Fill Restorative



Refill	Item #	Item #
Shade	Syringe	Capsule
A1	4866A1	4867A1
A2	4866A2	4867A2
A3	4866A3	4867A3
B1	4866B1	4867B1
C2	4866C2	4867C2

SYRINGE SHADE REFILLS:

Include 1 – 4g syringe; Technique Guide; Instructions for Use.

CAPSULE SHADE REFILLS:

Include 20 – 0.2g capsules; Technique Guide; Instructions for Use.

REFILLS—3M™ Filtek™ Supreme XTE Universal Restorative

Refill	Item #	Item #
Body Shades	Syringe	Capsule
A1B	4910A1B	4911A1B
A2B	4910A2B	4911A2B
A3B	4910A3B	4911A3B
A3.5B	4910A3.5B	4911A3.5B
B1B	4910B1B	4911B1B
B2B	4910B2B	4911B2B
C1B	4910C1B	4911C1B
C2B	4910C2B	4911C2B
A4B	4910A4B	4911A4B
A6B	4910A6B	4911A6B
B3B	4910B3B	4911B3B
B5B	4910B5B	4911B5B
C3B	4910C3B	4911C3B
D2B	4910D2B	4911D2B
D3B	4910D3B	4911D3B
White	4910WB	4911WB
Extra White	4910XWB	4911XWB

The 8
“Go-To”
Body
Shades

Dentin Shades	Syringe	Capsule
A1D	4910A1D	4911A1D
A2D	4910A2D	4911A2D
A3D	4910A3D	4911A3D
A4D	4910A4D	4911A4D
B3D	4910B3D	4911B3D
C4D	4910C4D	4911C4D
White	4910WD	4911WD

Enamel Shades	Syringe	Capsule
A1E	4910A1E	4911A1E
A2E	4910A2E	4911A2E
A3E	4910A3E	4911A3E
B1E	4910B1E	4911B1E
B2E	4910B2E	4911B2E
D2E	4910D2E	4911D2E
White	4910WE	4911WE
Extra White	4910XWE	4911XWE

Translucent Shades	Syringe	Capsule
Clear	4910CT	4911CT
Blue	4910BT	4911BT
Grey	4910GT	4911GT
Amber	4910AT	4911AT

SYRINGE SHADE REFILLS:

Include 1 – 3g syringe; Instructions For Use.

CAPSULE SHADE REFILLS:

Dentin, Enamel and Translucent Shade Refills include 10 – 0.2g capsules; Body Shade Refills include 20 – 0.2g capsules; Instructions for Use.

www.3M.com



3M Oral Care
2510 Conway Avenue
St. Paul, MN 55144-1000 USA
Phone 1-800-634-2249
Web 3M.com/dental

3M Deutschland GmbH
Location Seefeld
3M ESPE • ESPE Platz
82229 Seefeld • Germany
Info3MESPE@mmm.com
www.3MESPE.com

3M, Filtek and Z100 are trademarks of 3M or 3M Deutschland GmbH. Used under license in Canada. © 3M 2017. All rights reserved.