



TePe EasyPick™

Product Plastic core with cleaning surface of silicone.

Pack 12 pcs/blister pack, 36 pcs/blister pack, 48 pcs/blister pack, 60 pcs/blister pack. 10 blister packs/carton.

The clam pack (blister) is made of polyethylene terephthalate (PET). The backing card is made of cardboard and the folder of paper. The pocket case is made of polypropylene (PP).

The retail carton is made of chlorine free recyclable cardboard and the transport carton is made of corrugated cardboard.

Also available as 2-pack samples in plastic flow pack made of polypropylene (PP). Either 300 samples in a carton made of chlorine free recyclable corrugated cardboard, or 100 samples in a dispenser box made of FSC certified cardboard.

Recycling Product in household waste.

The blister pack should be separated and recycled. The plastic pack and pocket case in container for plastic waste and the backing card in container for recycling of cardboard.

The flowpack for samples in container for plastic.

The dispenser box in container for cardboard/corrugated cardboard.

The cartons in container for cardboard/corrugated cardboard.

Product features

- Easy cleaning between your teeth.
- 3 sizes to fit both narrow and wider interdental spaces.
- Ergonomic handle for easy use, non-slip grip.
- Efficient silicone cleaning surfaces.
- Durable.
- Developed in collaboration with dental experts.

